


PowerChek™ Pork Gelatin Kit Series



Halal Test in Gelatin Sample



Pork Gelatin Kit *from Extraction to Detection*

Reliable solutions for gelatin DNA extraction and pork detection


- Pork gelatin
- Marshmallow
- Gelatin processed food
- Gelatin raw material
- Various gelatin products

Halal Gelatin Testing

PowerChek™ Pork Gelatin Kit provides the state-of-the-art technologies for Halal Gelatin Testing from DNA extraction to detection, which can be checked to ensure its authenticity and that it is free from cross-contamination with porcine gelatin.




Work Flow








PowerChek™ Pork Gelatin Kit series

- Total solution kit ; From extraction to detection
- High sensitivity & specificity
- Internal control included
- Compatible with major Real-time PCR instruments



#	Sample
1	Pork raw meat
2	Pork noodle
3	Marshmallow (pork gelatin)
4	Porcine gelatin

Compatibility - Flexible open system

LT	BioRad	Roche	Agilent	Qiagen
ABI® 7500	CFX96™	LightCycler® 480	Mx3005P	Rotor-Gene Q®
				
●	●	●	●	●

ABI Prism® 7500 series are registered trademarks of ThermoFisher SCIENTIFIC. CFX96™ is a registered trademark of Bio-Rad Laboratories, Inc. Qiagen® is a registered trademark of Qiagen, Inc. ROCHE® is a trademark of Hoffman-LaRoche Inc. Rotor-Gene® is a registered trademark of Corbett Research Pty Ltd. Agilent Mx3000PTM and Mx3005PTM are a trademark of Agilent Technologies.

Ordering Information

Cat No.	Size	Product	Components	Applications
E0103	50	PowerPrep™ Gelatin DNA Extraction Kit	<ul style="list-style-type: none"> ▪ Lysis Buffer ▪ Elution Buffer ▪ Binding Buffer ▪ DNA Binding Column 	<ul style="list-style-type: none"> ▪ Proteinase K ▪ RNase A ▪ Wash Buffer
R0430E	50	PowerChek™ Pork Gelatin Real-time PCR Kit (FAM/HEX)	<ul style="list-style-type: none"> ▪ 2X Real-time PCR Master Mix ▪ Primer/Probe Mix ▪ Control DNA 	<ul style="list-style-type: none"> ▪ Gelatin, gelatin processed food ▪ Marshmallow, hard and soft capsules ▪ Various samples including gelatin ▪ Meat, meat products ▪ General processed food